

# How to Install TDO, Version 11 Best Practices

## How to Install TDO, Version 11 Best Practices

- Workstation should have Windows UNC set to low or turned off
- Domain Users needs to be a member of local administrators
- Browse to the TDO map network drive. Find folder “TDO Builds” and run the installer.
  - When running the installer, select “All defaults” and for the install type, select “Workstation.”
- After installer completes, it will create an icon on the desktop to launch the application. Right-click on the icon and select “Properties,” and make sure that it is set to run as administrator
- Display properties: monitor resolution needs to be 100% under the drop-down for "change the size of text, apps, and other items" or you can have problems with MS Word integration.
- <https://kb.tdo4endo.com/article.php?id=294>
- <https://kb.tdo4endo.com/article.php?id=107>
- <https://kb.tdo4endo.com/article.php?id=292>
- <https://kb.tdo4endo.com/article.php?id=340>
- Configure settings for that workstation under Administrative / Edit Control Table / Workstation Settings
- If the client uses CareStream, there is a DLL that TDO needs to update for the X-rays.

### Notes on monitor settings:

- The MS Word intergration for doing letters can have problems if the monitor is not set to 100% under display properties / "Change the size of text, apps, and other items"

## Notes on Windows Firewall Settings:

- If Windows firewall is enabled on the workstation, you have to create 2 firewall rules in the Windows Firewall. Both are outbound. One for TCP port 4150 and one for UDP port 4150

## DSLR Camera Notes - this is the camera attached to the microscope:

- You need to install CareStream driver SetupRVG-Carestream\_4.6.9.0-C.exe

Online URL: <https://kb.naturalnetworks.com/article.php?id=471>